



## **Nilotech Product Range**

Quality control instrumentation for specialized laboratories.



# Save Precious Time on Quality and Process Control

"When I use my NiloChecker, it is easy and fast to control gas and temperature in incubators. I simultaneously measure both parameters and quickly move to another incubator; the unit even saves all measurements. I always bring the NiloChecker into our various laboratories - it's easy to carry around and super accurate. A true all-in-one device."

## Kirsten Simonsen, MSc Laboratory manager, senior embryologist Maigaard Fertilitetsklinik

# Innovative, Exceptional and Dependable



"I am extremely satisfied with my NiloChecker tablet and sensor collection. In my experience, Nilotech is a passionate and engaged company that is committed to quality and their customers, and their Nilochecker product line clearly reflects this.

Despite being a high-performance tool, the Nilochecker is approachable, straightforward, and easy to use. It is robust and well-adapted to its lab application, and key design features help facilitate and ensure correct and accurate use. In cases where shortcomings have come up in specific end-use scenarios, Nilotech has welcomed the end-user feedback and worked to address them. They are continually refining and improving their designs and also coming up with new offerings. I am always excited to see these new items and put them to use in the field.

My own NiloChecker kit goes with me everywhere. Its accuracy, reliability, and utility are unparalleled and make it indispensable. Whether it's in my carry-on for the plane or in my pack for a local service visit, it's always in my bag and ready for use. Having used my Nilochecker extensively in the past few years, I can't imagine working without it now. It is an excellent tool that I depend on and use with confidence. Thank you, Nilotech!"

## Tim Carter, PhD Scientific and technical director Fertitech Canada

## Nilotech measuring probes











## NiloChecker 500

#### **One Reference instrument many applications**

The New all-in-one instrument made by Lab equipment specialists. NiloChecker guides you through control and calibration of critical equipment in your laboratory. No need for multiple instruments from various manufacturers. With the NiloChecker you can control and calibrate Incubators, IVF workstations, warming system etc. An increasing range of measurement probes ensures your investment in years to come.

### **Functions & Technology**

#### **Time saving**

Save precious time on process and quality control. Perform up to 10 simultaneous measurements. For instance, check your incubator for  $CO_{2^{\prime}}O_{2^{\prime}}$  and Temperature at the same time.

#### Easy as 1-2-3

A large touch screen with intuitive software guides you through the measurement. It helps you to perform measurements in accordance with international standards and manufacturers recommendations.

#### Integrated documentation

The Add-to-table function helps you processing and organizing data from many measurements.

#### Long battery life

Just bring the NiloChecker around. The high capacity battery ensures hours of measurement time between charging.

#### No need for calibration

The instrument is fully digital and will never need calibration. The instrument stays in the lab and measurement probes can be controlled and calibrated separately.

#### Compatible sensors

Compatible with  $CO_2$ ,  $O_2$ , temperature, and air velocity sensors. simultaneous measurements up to 10 sensors on 5 ports.

#### Integrated functions

Graph function: Important tool to track any measurement. Instant measuring: Instrument starts to measure once a probe is connected. Timer function: Start/Stop allows to measure in a predefined time. Calculations: Max, Min, Average, Moving average, Deviation from average. Easy to collect data from multiple sensors.

### Specifications

Measure range and accuracy Display and Touch Data storage and export Compliance Operating conditions Power supply Connectors Battery time Charging Time Depends on connected probe 7" - 800 x 480 pixel - Capacitive touch - Adjustable brightness Internal memory. - Export data to PC via USB ISO/EN 61010:2010 - IEC 60950-1:2005/ - AMD2:2013 - CE 0-50°C - 5-95% RH, IP 50 12V / 30W - Wall plug adaptor - 100-240V

5 Pcs - Stainless steel with lock and quick release

8 Hours - Level displayed in %

< 3 hours





Ordering info Part no: 111s001 - NiloChecker 500





## DG126 Digital CO<sub>2</sub>/O<sub>2</sub> probe with pump

DG126 Digital Carbon dioxide  $(CO_2)$  and Oxygen  $(O_2)$  probe with pump is revolutionary in design and performance. Due to the small design and wide connection possibilities, the probe is equally useful for controlling box-type and desktop incubators.

## **Functions & Technology**

#### **Dual measurements**

Connect the flexible tubing to a sample port or use the sample needle to sample from narrow openings. Start the pump from the NiloChecker device, and it measures  $CO_2$  and  $O_2$  simultaneously.

#### Flexible design

The probe can be placed near- or on- the target device to minimize tubing length. Short tubing minimizes measurement time and usage of gas.

#### Ease of use

Connect up to 5 DG126 probes simultaneously. With the NiloChecker in your hand, you can start/stop measurements, monitor progress and store readings for later analysis.

#### Output data

 $\rm CO_2$  and  $\rm O_2$  in %. The NiloChecker calculates and maps data including min., average, and max. values from several measurement sessions. See NiloChecker product sheet for more information about exporting measurement data.



## **Specifications**

Sensors Measure range Accuracy T90 Compensation for ambient conditions Pump flow Display resolution and update Cable lenght Tube fitting Compliances Calibration NDIR (CO<sub>2</sub>), Electro-chemical (O<sub>2</sub>) 0-30% O<sub>2</sub> / 0-12% CO<sub>2</sub> 0,2% O<sub>2</sub> , 0,1% CO<sub>2</sub> +/- 3% of reading <20 sec. (O<sub>2</sub>), < 15sec (CO<sub>2</sub>) Temperature (20-40°C), altitude (700-1100 mbar) 100-150 ml/min 0,1% ,1 update per sec. 2 m Luer-lock as standard. Other 1/8" fittings optional. CE with NiloChecker 500. RoHs. Delivered with factory calibration certificate.. Can be calibrated in accordance with ISO/IEC 17025. 0-50°C, Calibration requires adaptor. Part no: 115s001. 5-95% RH (Non condensing)





Operating conditions



Ordering info

 Part no: 126s001 - DG126 Digital CO2/O2 probe with pump.

 Delivered with factory calibration certificate.

 Part no: 800s004 - Calibration of Nilotech gasprobe incl. factory and gas certificates

#### Nilotech Aps.





## DG112 Digital CO<sub>2</sub>/O<sub>2</sub> probe

Measuring gas inside incubators are tricky, especially with gas analyzers aspirating gas from the incubator. Incubator stability can be influenced by even minute amounts of gas aspiration and the question often arise: When has enough gas been sampled. Measuring inside the incubator instead of aspirating gas, often result in more precise measurements.

## **Functions & Technology**

CO<sub>2</sub>/O<sub>2</sub> probe

The new Digital Carbon dioxide (CO<sub>2</sub>) and Oxygen (0<sub>2</sub>) probe is revolutionary in design and performance.

Due to the flat design and flexible connection, the probe can be inserted directly into the desktop incubator chamber, allowing for long-term continuous measurements and not only spot-check. See NiloChecker documentation for advanced functions





## **Specifications**

Sensors Measure range Accuracy T90 Compensation for ambient conditions Display resolution and update Cable lenght Compliances Calibration

NDIR (CO<sub>2</sub>), Electro-chemical (O<sub>2</sub>) 0-30% O<sub>2</sub> / 0-12% CO<sub>2</sub> 0,2% 0, 0,1% CO, +/- 3% of reading <20 sec. (O<sub>2</sub>), < 15sec (CO<sub>2</sub>) Temperature (20-40°C), altitude (700-1100 mbar) 0,1% ,1 update per sec. 2 m CE with NiloChecker 500. RoHs. Delivered with factory calibration certificate.. Can be calibrated in accordance with ISO/IEC 17025. 0-50°C, Calibration requires adaptor. Part no: 115s001. 5-95% RH (Non condensing)

Operating conditions

Part no: 112s003 - DG112 VFC (Very flat cable for desktop incubators) Delivered with Factory calibration certificate. Part no: 800s004 - Calibration of Nilotech gasprobe incl. factory and gas certificates







Nilotech Aps.

Knudstrupvej 14 - DK-4270 Hoeng, Denmark +45 30 32 32 96 - contact@nitotech.eu - www.nilotech.eu Document 112d002-2





## DT112 Digital temperature probe

Correct temperature inside incubators or on heating plates is crucial for the survival of living cells. The flat temperature probe can be inserted into any type of incubator, without disturbing conditions. The cobber bottom ensures good thermal contact to any type of heated surface and including non-metal surfaces.

## **Functions & Technology**

#### Temperture probe

This new surface-temperature probe is ideal for controlling and calibration incubators and heating surfaces.

The bottom is made from copper to give optimal thermal contact. The probe can be placed on incubator shelf or bottom. Or it can be placed inside dishware. The insulated top minimize disturbance, when used for heating plates or inside flow cabinets.





## Specifications

Sensor type Measure range & accuracy Display resolution and update T90 (min:sec) Compliance

> Operating conditions Cable lenght Materials

> > IP Class Calibration

-20-50°C / 0,03°C accuracy possible. 0,01 °C, 1 update per second 1:45 ISO/EN61010-1 - RoHS Medical grade sensor - ISO 80601-2-56 0-50°C / 5-95% RH 1,7m (including 0,5m flat cable) Housing: ABS, Copper. Cable: TPE (Flat part), PVC (Round part) IP50 Delivered with factory calibration certificate. Accuracy +/- 0,1°C Can be calibrated in accordance with ISO/IEC 17025. Calibration requires adaptor. Part no: 115s001.

Digital temperature sensor





Ordering info

Part no: 112s001 - DT 112 FC Digital Temperature Probe Part no: 800s002 - Accredited calibration of 1 temperature sensor

> Nilotech Aps. Knudstrupvej 14 - DK-4270 Hoeng, Denmark +45 30 32 32 96 - contact@nitotech.eu - www.nilotech.eu

Document 112d001





## DT122 Dual digital temperature probe

Correct temperature inside incubators or on heating plates is crucial for the survival of living cells.

DT122 Dual Digital Temperature Probe is ideal for controlling and calibrating desktop incubators, box-type incubators and heated surfaces.

## **Functions & Technology**

#### Dual temperature probe

DT122 has 2 separate temperature sensors in one probe. With NiloChecker 500, you can measure and calibrate up to 10 positions or chambers simultaneously.

The unique sensor design has a copper bottom and insulating top, which ensures optimal thermal contact. All our sensors comes factory calibrated and can be delivered with a 3rd party certificate of accredited calibration.





## Specifications

Sensor type Measure range & accuracy Display resolution and update T90 (min:sec) Compliance

> Operating conditions Cable lenght Materials

> > IP Class Calibration

2 x Digital temperature sensors -20-50°C / 0,03°C accuracy possible. 0,01 °C, 1 update per second 1:45 ISO/EN61010-1 - RoHS Medical grade sensor - ISO 80601-2-56 0-50°C / 5-95% RH 1,7m (including 0,5m flat cable) Housing: ABS, Copper. Cable: TPE (Flat part), PVC (Round part) IP50 Delivered with factory calibration certificate. Accuracy +/- 0,1°C Can be calibrated in accordance with ISO/IEC 17025. Calibration requires adaptor. Part no: 115s001.







 Part no: 122s001 - DT122 Dual Digital Temperature Probe

 (Delivered with factory calibration certificate).

 Part no: 800s003 - Accredited calibration of 2 temperature sensors.

#### Nilotech Aps.

Knudstrupvej 14 - DK-4270 Hoeng, Denmark +45 30 32 32 96 - contact@nitotech.eu - www.nilotech.eu Document 122d001-1





## DT121

## Up-Down digital temperature probe

Ensuring the correct temperature inside incubators is crucial for the survival of living cells. Our new DT121 Up-Down Digital Temperature Probe has a unique design with a copper bottom and top which is separated by a spring. This ensures the perfect thermal contact to chamber floor and lid. For the first time you can measure both the top and chamber floor temperature at the same time with only one probe.

## **Functions & Technology**

#### **Dual measurements**

Our unique Up-Down probe is ideal for controlling and calibrating desktop incubators through its 2 digital temperture sensors.

#### Flexible design

The probe can be placed on a incubator chamber floor, while the spring ensures good contact between upper sensor and lid.

#### Ease of use

We have made it possible for you to connect up to 5 probes in one NiloChecker. Controlling and calibrating desktop incubators has never been faster og easier.



## Specifications

#### Sensor type

Measure range & accuracy Display resolution and update T90 (min:sec) Compliance

> Operating conditions Cable lenght Materials

> > IP Class Calibration

2 x Digital temperature sensors -20-50°C / 0,03°C accuracy possible. 0,01 °C, 1 update per second 1:45 ISO/EN61010-1 - RoHS Medical grade sensor - ISO 80601-2-56 0-50°C / 5-95% RH 1,7m (including 0,5m flat cable) Housing: Copper. Cable: TPE (Flat part), PVC (Round part) IP50 Delivered with factory calibration certificate. Accuracy +/- 0,1°C Can be calibrated in accordance with ISO/IEC 17025. Calibration requires adaptor. Part no: 115s001.





Ordering info

Part no: 112s001 - DT 121 Up-Down Digital Temperature Probe (Delivered with factory calibration certificate) Part no: 800s003 - Accredited calibration of 2 Temperature sensors.

#### Nilotech Aps.

Knudstrupvej 14 - DK-4270 Hoeng, Denmark +45 30 32 32 96 - contact@nitotech.eu - www.nilotech.eu Document 112d001\_2



## DT123 Digital dish-temperature probe

Ensuring the correct temperature inside dishes is crucial for the survival of living cells. Our new DT 123 Digital dish-temperature probe has a very thin 200 mm flexible wire and a small sensor at the tip. (4 mm wide, 0,5 mm thickness) This design allows you to measure liquid temperature in dishes or tubes. Due to the slim design the sensor temperature stabilizes quickly.

## **Functions & Technology**

## Dedicated temperature sensor for dishes and tubes.

DT 123 is ideal for controlling and calibrating heated stages and other sensitive heating devices. The probe can be placed on any surface, such as metal, glass, plastic etc. We have made it possible to connect 5 different probes to one NiloChecker. Controlling and calibrating your lab equipment has never been faster or easier.





## **Specifications**

Sensor type Measure range & accuracy Display resolution and update T90 (min:sec) Compliance

> Operating conditions Cable lenght Materials

> > IP Class Calibration

Digital temperature sensor -20-50°C / 0,03°C accuracy possible. 0,01 °C, 1 update per second 0:45 ISO/EN61010-1 - RoHS Medical grade sensor - ISO 80601-2-56 0-50°C / 5-95% RH 1,7m (including 0,5m flat cable) Flexible FR-4, Epoxy. Cable: TPE (Flat part), PVC (Round part) IP50 Delivered with factory calibration certificate. Accuracy +/- 0,1°C Can be calibrated in accordance with ISO/IEC 17025. Calibration requires adaptor. Part no: 115s001.



Ordering info

 Part no: 123s001 - DT 123 Digital dish temperature Probe

 (Delivered with factory calibration certificate)

 Part no: 800s002 - Accredited calibration of 1 temperature sensor.

#### Nilotech Aps.

Knudstrupvej 14 - DK-4270 Hoeng, Denmark +45 30 32 32 96 - contact@nitotech.eu - www.nilotech.eu Document 123d001\_2







### DI V119

## Digital low velocity air flow probe

Correct laminar flow in IVF workstations is crucial to avoid contaminating cells and ensuring operator safety. The DLV119 can be used to control filters and flow according to manufacturer's recommendations and international standards such as EN 12469 in:

- Biological Safety Cabinets Chemical Fume Hoods
- Laminar Flow Hoods
- Clean Benches
- HEPA & Filter Boxes
- And other

DLV119 is a versatile and rugged, high-performance air velocity probe for low flow.

## **Functions & Technology**

#### **Functions and Technology**

DLV119 is suitable for demanding applications, including those in corrosive or alkaline environments. With its robust, splash proof design, the DLV119 is designed to support a wide range of product and process control air flow applications.

Software in the NiloChecker supports calculation of airflow in accordance with international standards and easy reporting.





## Specifications

Sensors Measure range Accuracy Response time Compensation for ambient conditions Display resolution and update Cable lenght Compliances Calibration

> Operating conditions Materials

Flow: 0,15 - 1,5 m/s. Temperature:  $0^{\circ}$ C -  $60^{\circ}$ C Flow:  $\pm 1\%$  of reading + 0,05 m/s. Temperature: +/-1°C 400 ms Temperature (0-60°C) 0,01 m/s, 0,01°C,1 update per sec. 2 m CE with NiloChecker 500. RoHs. Delivered with factory calibration certificate.. Can be calibrated in accordance with ISO/IEC 17025. Calibration requires adaptor. Part no: 115s001.

Hot wire anemometer. Temperature compensated.

5-95% RH (Non condensing) Cable: PVC coated

Housing: Polycarbonate (PC), UL94-V0 (head) UL94-HB (housing). Aluminum (Cable ring)





Ordering info

Part no: 119s001 - DLV119 Digital Air Velocity Probe (Delivered with factory calibration)



Document 119d001





## **NiloBlock**

#### A new heating block concept for oocyte retrieval

A new insulated heating block which brings optimal temperature stability to test tubes. This newly developed concept allows visual inspection of tube content with minimal heat loss. NiloBlock can be used with standard heating surfaces and compatible block heaters.

## **Functions & Technology**



Materials

Weight Contry of origin

Cleaning

Block: Anodized aluminum Cover: Polycarbonate 630g Denmark

Block: Handwash and/or autoclave Cover: Handwash or machinewash (max 70°C)





Ordering info Part no: 110S002 - NiloBlock



## Accessories

Accessories for our product range



**Calibration adapter** Adaptor to check and adjust any NiloChecker probe.

**Ordering info** Part no: 115s001



Water trap for DG126 Water trap with luer-lock fittings.

**Ordering info** Part no:



**Calibration cap for DG112** Silicone cap for callibrating the DG112 - CO2 /O2 probe.

**Ordering info** Part no: 112s002



**Extension cable** 3 m. Extension cable for NiloChecker probes

**Ordering info** Part no: 112s005



NiloCase Mobile hard case for NiloChecker and sensors.

**Ordering info** Part no: 111s020







## Accessories

Accessories for our product range

## **UMONIUM38® NEUTRALIS TISSUES**

NEUTRALIS TISSUES are broad-spectrum cleaning and disinfectant wipes. They can be used to disinfect your NiloChecker, NiloChecker probes and other non-invasive medical devices. This unperfumed and uncolored product has been developed for use in neonatal incubators, IVF laboratories and more. Delivered in container with 100 wipes 20x20 cm



**Ordering info** Part no: 130s001

